NEATH PORT TALBOT COLLEGE COLEG CASTELL NEDD PORT TALBOT

School of Maths & Science Science Practical

Blood Typing

♦ Aim

At the end of this experiment you should be able to:

- 1. Examine different blood types and show how they can be used to work out an unknown blood type.
- 2. Analyse the results and assess their value as evidence in court.

♦ Introduction

ABO Blood Typing

Once it has been established that the stain at a scene is blood (Kastel-Meyer) and that it is human (Precipitin), it can then be further established whether the blood relates to a particular person. This test relies on the reaction of antisera with antigens. There are four main blood groups, A, B, O and AB based on the presence or absence of antigens on the surface of the blood cells.

♦ Safety

Control Measures

- The wearing of safety **glasses** and **a laboratory coat at** all **times** will be sufficient to take account of most hazards and significant risks
- All waste is to be placed in the labeled container immediately after use
- You are reminded of the need of good laboratory practice in order to maintain a safe working environment.

Hazards

♦ Procedure

- 1. Collect four plastic drop plates, three different coloured mixing sticks and an unknown blood sample.
- 2. Using a pipette, drop 3 5 drops of the unknown blood onto each of the three wells in one of the drop plates.
- 3. To the well marked A, add the Anti-A serum, to the well marked B add the Anti-B serum and to the well marked Rh add the serum marked Anti-Rh.
- 4. Use a different coloured stick to mix the serums into the blood drops for a few minutes.
- 5. The blood may change consistency, clumping together and turning opaque, this is known as agglutination and is positive for the serum used.
- 6. If agglutination occurs in the well marked A then it is positive for type A blood, If it occurs in B then it is positive for type B, if it agglutinates in both A and B then it is positive for type AB and it is considered to be type O if neither A nor B agglutinates.
- 7. If agglutination occurs in the Rh well then it is Rhesus positive (+), no change indicates rhesus negative (-).
- 8. Once the type of unknown blood has been established you can compare it to the three suspect or elimination samples provided.
- 9. For each of the blood types, follow the process outlined above until a match is made.

Sample	Anti A	Anti B	Rhesus	Result
1				
2				
3				
Unknown				

♦ Results

